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When exposed to a long use and high temperatures, frying oils used on human feed, decomposed, and their properties become altered by the subjection to these elements.

As long as an oil is used to deep fry some food products, it increases the amount of Total Polar Compounds (TPC) which can determinate the appearance of health nasty problems such as cancer, all along the human digestive system.

Some of the products obtained as the result of the continuous oil heating, have mutagenic - carcinogen - properties.

So, it's fundamental to implement quality control procedures in all the places where fry food is made, in order to assess, at each moment, the quality level and status of the oil.

The degradation level can be assessed by an **easy-to-use, quick, accurate and reliable** method - - conceived to be used by all food operators - consisting on a three steps color change system, with the possibility to an immediate reading and comprehension of the test result.

According to the concentration level of the Total Polar Compounds in the oil, the color changes from blue to olive green and all the operator needs to do is to compare the obtained color of the product with the scale on the plastic tube.

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

**INSTRUCTIONS**

**1** - Heat the oil around 120°C (high enough to melt the product inside the tube). Fill the tube 'till the lower label limit, using the filler (supplied with the kit).

**2** - Cap the tube and shake well ('till homogenize the mixture).  
**Note:** If the oil is not warm enough it's possible to heat the tube in a warm water bath ( $\pm 65^\circ\text{C}$ ) avoiding to wet the label.

**3** - Wait 2 minutes and hold the tube on the counter-light. Compare with the chromatic scale on the label.

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**ASSESSMENT SCHEDULE**

QUALITY (according to the Results)	Esteemed Concentration of TPC	Correct procedure
GOOD	0 - 5%	USE WITH RELIABILITY
MEDIUM (close to GOOD)	6 - 11%	USE WITH RELIABILITY
MEDIUM (close to BAD)	12 - 17%	USE WITH RELIABILITY (*)
BAD	> 17%	REJECT

(\*) This oil must be controlled after the next fry process

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Kit must be stored in a dry and cold place, preserved from sunlight.